



What's New Questionnaire August 2018

Do you have questions on the results of our WNQ? If so, please reach out to our ABO certified Customer Service team (920) 469-2559. Additionally, if you have any suggestions for topics on next month's WNQ - please don't hesitate to let us know! Contact sara@cherryopticalinc.com today with your best trivia questions!

"What you do today can improve all of your tomorrows." - Ralph Marston

<3 the Cherry Optical, Inc Team

Congratulations to:

1) True or False? When ordering any of the Eyezen+ lenses an OC height is required.

2) True or False? When measuring the Eyezen+ 0, it's measured at the 180* line.

3) A +3.00D lens has the following prismatic effect 8 mm below the optical center:

- A.) 2.4 base in
- B.) 2.4 base out
- C.) 2.4 base down
- D.) 2.4 base up

4) True or False? When measuring the Eyezen+ 1, 2 or 3 the power is measured at 4 above.

5) If the pads are brought closer together, this adjustment will:

- A.) Drop the front
- B.) Raise the front
- C.) Set the front closer to the face
- D.) Loosen the temple tension

6) When light entering the eye, without corrective lenses, comes to a focus in front of the retina, the patients refractive error is:

- A.) Myopia
- B.) Presbyopia
- C.) Hyperopia
- D.) Antimetropia